

15000 SHELDON SCIENCE/MATH DISCOVERY CENTER

47-1/4" L. X 47-1/4" W. X 26-1/2" to 36-1/2" H. in 2" increments. The imprints on the top in English and Metric provide the mathematical bases for math and for quantitative treatment of science. Vertical scales provide the third dimension. Center includes imprinted work top with square foot/square inch grids; 3-decimeter/centimeter grids; English/Metric scales; 360 degree protractor with clockwise/counter-clockwise scales, and extended protractor points. Write on top of tan plastic, suitable for soap base crayons and similar markers. Raised edge on top to keep equipment from rolling off and facilitate cleaning-up. Unit furnished with two (2) complete sets of easy on/off C-clamp rods, two (2) carryalls and two (2) totetrays. Adjustable metal legs to be 14 gauge with stainless foot cover.