

15100 SHELDON MOBILE STUDENT SCIENCE LEARNING CENTER

48" L. X 24" W. X 25" H. Drop leafs are 15" wide. Unit is constructed of 7-ply hardwood plywood. Top is 3/4" imprinted grid plastic laminate with black "T" molding all exposed sides. Center furnished complete with the following: two (2) easy/on, easy/off c-clamp rods, one (1) crossbar, two (2) greenlaw arms, two (2) connectors, eleven (11) white fiberglass tote trays (5" X 3" X 19" overall), three (3) fiberglass tote trays, (17" X 11-1/2" X 3-3/4") four (4) locking swivel castors (shipped separate and field installed). Center comes complete with locks.