

A. **Sheldon 20625 Instructor Demonstration Center:**
60" L. X 31" W. X 36" H. Sink on Right End.

TP 12.14.17

1. Demonstration Center: Includes one (1) S183528-100 sink cabinet, one (1) B183528-571 five-drawer cabinet, one (1) kneespace, and a finished back. The Center also includes a black Shelresin top with epoxy resin drop-in sink on right end. **Internal service piping and wiring not included.**
2. Standard services and accessories included:
 - a. One (1) 80020-CV Unicast combination cold water and double gas fixture with check valve.
One (1) D05 Shelresin 14" x 10" x 6" I.D. drop-in sink.
One (1) sink outlet, strainer and stopper.
One (1) 85101 GFI duplex electrical outlet.
One (1) 86380 upright rod assembly with two (2) 86320 rod bases.