

A. **Sheldon 20655 Instructor Demonstration Center:**  
96" L. X 31" W. X 36" H. Sink on Right End.

TP 12.14.17

1. Demonstration Center: Includes one (1) S353528-200 sink cabinet, one (1) B353528-570 seven-drawer cabinet, one (1) 24" kneespace, and a finished back. The Center also includes a black Shelresin top with epoxy resin drop-in sink on right end. **Internal service piping and wiring not included.**
2. Standard services and accessories included:
  - a. One (1) 80020-CV Unicast combination cold water and double gas fixture with check valves.  
One (1) D20 Shelresin 16" x 16" x 7-1/2" I.D. drop-in sink.  
One (1) sink outlet, strainer and stopper.  
One (1) 85101 GFI duplex electrical outlet.  
One (1) 86380 upright rod assembly with two (2) 86320 rod bases.