

## **20692 PORTABLE DEMONSTRATION CENTER**

49"L. x 24"W. x 36"H. Unit features a (1") thick Phenolic resin countertop; two (2) hinged doors. Assembly is provided with finished back. Unit includes: four (4) 4" rubber tired swivel casters, two (2) with locking brakes. Service includes a 1-gallon polyethylene water supply container connected to a chrome-plated galley pump. Stainless steel sink connected to a 1-gallon polyethylene waste container. Four (4) plastic tote trays (16 ½" x 12" x 4"D). One (1) 110V G.F.I. duplex receptacle with 10'L., 3-wire extension cord and grounded plug. One upright rod assembly set.

Rev. 11.23.2011