

20693 SHELDON MOBILE DEMONSTRATION TABLE

49" L. X 24" W. X 36" H. Mobile Instructor's Demo Center features a 49" L. X 24" X 1" Phenolic Resin Top. Mobile cabinet with two hinged doors, moves freely on four (4) swivel casters, two (2) with locking brakes. Services include 1-gallon polyethylene water supply container connected to chrome-plated rocker galley pump. Stainless steel sink is connected to 1-gallon polyethylene waste container. Containers are easily removable. One duplex electrical receptacle with 10' extension cord and grounded plug, four plastic tote trays (16½" x 12" x 4"), one (1) #86320 upright rod socket and one (1) #86385 upright rod assembly.

Rev. 10.12.11