

A. **Sheldon 20698 ADA Instructor Demonstration Center:**  
96" l. x 31" w. X 33" H. Sink on Right End.

TP 12.14.17

1. Demonstration Center: Includes one (1) S343228-000 ADA sink cabinet with removable panel, one (1) B303228-560 six-drawer cabinet, one (1) ADA kneespace, and a finished back. The Center also includes a black Shelresin top with epoxy resin ADA drop-in sink on right end. **Internal service piping and wiring not included.**
2. Standard services and accessories included:
  - a. One (1) 80020-WB-CV Unicast combination cw/gas fixture with wrist blade handle and check valves.  
One (1) A25 Shelresin 18" x 15" x 5" I.D. drop-in sink.  
One (1) sink outlet, strainer and stopper.  
One (1) 85101 GFI duplex electrical outlet.  
One (1) 86380 upright rod assembly with two (2) 86320 rod bases.