

20795 SHELDON INSTRUCTOR DEMONSTRATION CENTER

68" L. X 24" W. X 36" H. Mobile Instructor's Demo Unit provided with one (1) storage cubical (15"D. x 20"L. x 16-3/4"H) with three drawers. Top is 1" Shelresin with one edge angled to interface with the Sheldon Trifacial Service Island or wall sink. Assembly is provided with finished back; four (4) rubber-tired swivel casters, two (2) with locking brake; 110V electrical raceway with three plug-in receptacles; full length modesty panel; and one (1) 86375 clamp on upright rod assembly. Rod storage is incorporated into the frame.

Rev. 10.24.11