

## **27860 SHELDON RINSEAWAY SINK STATION (Left Sink)**

70-1/2"L. X 24"W. X 29"H. to counter top. Consisting of one piece sloped counter and integral sink (minimum 18" X 12" X 5" I.D.) located on the left with raised edges on all sides and coved corners; top to be of fiberglass reinforced structural polyester.

Services include one (1) 80022 Unicast cold water fixture, one (1) G-5020 hand-held drench hose, one (1) 35" base cabinet with three (3) half-width drawers, cupboard section with one (1) hinged door and one (1) adjustable shelf. One (1) 35" sink cabinet with two (2) hinged doors.

Rev. 11.23.2011