

**27865 SHELDON RINSEAWAY SINK STATION**

70-1/2" L. X 24" W. X 36" H. to counter top. Consisting of one-piece sloped counter and integral sink (minimum 18" X 12" X 5" I.D.) with raised edges on all sides and covered corners; top to be of glass fiber reinforced structural polyester; services include cold water gooseneck fitting; one (1) base cabinet with four (4) halfwidth drawers, cupboard section with one (1) hinged door and one (1) adjustable shelf. One (1) sink cabinet with two (2) hinged doors.