

27870 SHELDON RINSEAWAY SINK STATION (Right Sink)

70-1/2"L. X 24"W. X 29"H. to counter top. Consisting of one piece sloped counter and integral sink (minimum 18" X 12" X 5" I.D.) located on the right with raised edges on all sides and coved corners. Top to be of fiberglass reinforced structural polyester.

Services include one (1) 80022 Unicast cold water fixture, one (1) G-5020 hand-held drench hose, one (1) 35" base cabinet with three (3) half-width drawers, cupboard section with one (1) hinged door and one (1) adjustable shelf. One (1) 35" sink cabinet with two (2) hinged doors.

Rev. 11.23.2011