

66100 SHELDON ADA SAFETY CENTER

35" L. X 22" D. X 85" H. Overall. Unit shall be constructed as detailed within "Casework Construction" section of this specification. Unit provided with two (2) gentle spray eyewash outlets located 34" above finished floor for individual in wheelchair accommodation. Eyewash assembly located over special depth and size molded fiberglass top/sink combination with raised edges all sides. Overhead shower modified to have proper distance from case to accommodate emergency individual and with proper pull handle a maximum of 48" above finished floor. Unit shall have one (1) fixed shelf for storage of fire blanket or first aid kit (not included). Unit shall have removable back panel, and a lower removable epoxy coated steel plumbing enclosure, properly contoured for required clearance to accommodate individual in wheelchair. Unit shall include a sign indicating the "International Symbol of Accessibility".

Rev. 1.31.12