

76810 SHELDON WALL SINK CENTER

Size 30"W x 74"L x 61"H Overall. Wall sink assembly shall consist of two supporting base cabinets, stainless steel tub sink with integral service turret and service fixtures, and service wall assembly.

Each base cabinet shall be a special 24"L unit with one end panel 29"H, and opposite panel 34-3/4" H. Cabinets 22" deep with 7" deep pipe chase with end closure panels. Each cabinet with a single hinged door and a full length fascia panel connecting both units. A removable closure panel shall be secured between cabinets allowing access to all plumbing services. Cabinet and closure panel shall be constructed identical to other units provided within this specification.

Stainless steel tub sink shall be Type 302 with integral 6"H back panel having an inside dimension no less than 72" x 22" x 8" deep. Entire tub sink shall be formed to drain to center outlet. Integral back panel provided with a formed stainless steel removable cap full length 7-1/2" wide. Wall clips provided for securing removable cap. Integral back panel provided with three panel mounted hot and cold mixing goosenecks with aerators.

Service wall assembly shall consist of a 3/4" thick, 72"L x 21"H plastic laminate back panel with an 8" deep, full length laminate shelf. Provide four panel mount K & V #1132 9" long utility hooks.