

76840 SHELDON ISLAND SINK STATION

49"L. x 48"W. x 36"H. to counter top. 1" plastic laminate top, furnished with two (2) 19 ½"L. x 19"W. x 7 ½"H. stainless steel sink with outlet, strainer, and stopper (bowl size 16" X 13 1/2" X 7 1/2 I.D.) Also furnished with two (2) 80100 combination hot and cold water gooseneck fixture. Two (2) sink cabinets, back to back, are 47" wide with 16" deep cupboards and (2) hinged doors.

Rev. 1.26.12