

SC7131 FLAMMABLE STORAGE CABINET

43"W. x 44 5/8"H. x 18"D. Cabinet is constructed of one-inch (1") thick, high-density, 9-ply, exterior grade plywood finished with yellow epoxy paint (inside and out) and 10" x 10" hazard labels on the doors. Unit furnished with continuous piano-type hinges on doors, lock with two keys, two (2) adjustable shelves mounted on metal shelf standards and supports, and corrosion-resistant hardware. Cabinet floor constructed as a two-inch (2") liquid-tight trough to contain spillage. Cabinet complies with all O.S.H.A. and National Fire Protection Association standards. Capacity: 45 Gallons

Rev. 1.31.12