

## **SC8051 WOOD ACID STORAGE CABINET**

31"W. x 36 5/8"H. x 20"D. Cabinet is constructed of one-inch (1") thick, high-density, 9-ply, exterior grade plywood finished with multiple coats of epoxy paint. Cabinet bottom constructed as a liquid-tight, two-inch (2") trough to contain accidental spills. Top is recessed to include a liquid-tight polypropylene tray. Tray will hold 3 gallons of spilled liquid and can be used as a dispensing area. Interior is fully lined with 1/8" thick polypropylene and all screws and fasteners are to be plastic. Cabinet has an interlocking door assembly, lock hasps and padlock with two keys' and one fixed shelf. The word "ACID" is printed on the cabinet in six-inch (6") red letters. Unit furnished with wooden door handle and four leveling feet. Cabinet complies with all O.S.H.A. and National Fire Protection Association standards. Capacity: (30) 2.5 liter or (90) 500ml bottles.

Rev. 1.31.12