

**SHELDON TWO-STUDENT TABLE**  
**U542130-I, U542430-I, U542136-I, U542436-I**

(54" X 21" X 30") - U542130-I

(54" X 24" X 30") - U542430-I

(54" X 21" X 36") - U542136-I

(54" X 24" X 36") - U542436-I

Overall - Hardwood fully framed with rails 3/4" X 4" with diagonal corner braces tenoned and screwed to inner face of rails. The corner braces shall be metal, grooved and screwed with four (4) screws. Table legs shall be properly fitted into position and securely fastened to diagonal corner brace with nut, washer and 3-1/2" x 5/16 carriage bolt completely running through the leg providing a positive system whereby bolt can be tightened without depending upon screw holding power of the table legs. Legs shall be 2-1/4" square of solid laminated hardwood, thoroughly glued, and corners rounded. Legs shall be equipped with rubber leg shoes, and nylon adjustable glides.

1" Imprinted Top - measuring units are permanently printed on the plastic sheet which is overlaid with clear melamine. Check points and computations can be made directly on the surface with a crayon or grease pencil, and a damp cloth wipes everything clean. Included are linear scales, grids and squares, diagonals, circles, protractors -- practically everything your students need for measuring, making graphs, charts and recording results.