

31125 SHELDON PLANT GROWTH CENTER

Overall 26-1/2" L., 30" W., 60" H., furnished with four (4) rubber tired casters, each with locking device. The single light fixture is complete with eight (8) fluorescent Grow-lux tubes with toggle switch control for each set of four (4) lights. A twenty-four hour timer is affixed into the wood upright post close to the light fixture, and has a pigtail connection cord. The metal enclosure stores the ballast and wiring for the lights. The light fixture is securely mounted to the posts. The storage shelves are 24" X 30". The top shelf is adjustable and sits 12" below the light tubes and can be lowered as plants grow to 18" below the lights. Unit contains six (6) fiberglass growing trays.