

31435 SHELDON ANIMAL STUDY CENTER

56" L. X 27" W. X 74" H. Finished oak hardwood exterior. Top superstructure of steel with baked synthetic finish. Top removable for access to wiring. Interior steel top and ends finished in white. Ends contain integral shelf supports on (1-1/2") centers and perforated to assure even air circulation regardless of shelf location. Water-tight stainless steel pullout debris trays at bottom are easily removed for cleaning from the front without disturbing interior conditions. Removable mesh bottom and three (3) adjustable plated wire shelves. Two (2) cupboard storage compartments for food, bedding and one vented compartment (for cooling unit, which is available when specified at an adjusted price) is in base. Two (2) sliding glass doors and instrument panel cover provided with locks. Switches on instrument panel control master circuit breaker switch, temperature from ambient up to 38 degrees Celsius (+ 1 degree ambient based on 22 degrees Celsius; down to 8 degrees Celsius if cooling unit is included), lighting and blower. Hinged instrument panel cover is transparent Plexiglas with lock. Forty (40) watt fluorescent lamp housed in top. Positive ventilation is insured by an integral blower and (3") diameter duct stub to vent outdoors. Blower to provide four air changes per minute. Contract rating 1,100 watts total heater load. Unit vents into room and includes specially formulated activated charcoal filter and dust filter, both easily accessible for servicing. Ducting to outside not required.