

## 31440 SHELDON ANIMAL STUDY CENTER

56" L. X 27" W. X 74" H. Finished oak hardwood exterior. Top superstructure of steel with black powder-coat epoxy finish. Top removable for access to wiring. Interior steel top and ends finished in white. Ends contain integral shelf supports on (1/2") centers and perforated to assure even air circulation regardless of shelf location. Water-tight stainless steel pullout debris trays at bottom are easily removed for cleaning from the front without disturbing interior conditions. Removable mesh bottom and three (3) adjustable plated wire shelves. Two (2) cupboard storage compartments for food, bedding and one vented compartment (for cooling unit, available when specified at an adjustable price) are in base. Two (2) sliding glass doors and instrument panel cover provided with locks. Switches on instrument panel control master circuit breaker switch, lighting, blower, and temperature (temperature ranges from ambient, based on 22 degrees Celsius, up to 38 degrees Celsius +/- 1 degree). Hinged instrument panel cover is transparent Plexiglas with lock. Thirty-two (32) watt fluorescent lamp housed in top. Positive ventilation is insured by an integral blower and (3") diameter duct stub to vent outdoors. Blower to provide four air changes per minute. Contract rating 1,100 watts total heater load. Unit vents into room and includes specially formulated activated charcoal filter and dust filter, both easily accessible for servicing. Ducting to outside not required. Unit mounted on heavy-duty casters and furnished with a finished tack-board back.

Rev. 10.5.11