

31480 EARTH SCIENCE STUDY CENTER

61"L. X 24"W. X 42"H. Center of formed plastic with all corners covered and with lower end forming the lake area for water depths up to (2 ½"). Clear acrylic contour map grid (24" X 25 ½") and marked with (2") squares. Plastic laminate form fits into the lake area which enables flooding conditions and permits the drawing and studying of contour maps. Center shall include sedimentation chambers (48" X 12" X 1 ½") with an inclined trough, crustal stress simulator (35 ½" X 13 ¼" X 3"), and rainfall simulator. Countertop is (1") Shelresin. Needs electrical source.

Rev. 10.5.11