



AXIS INFINITY®

The **Axis Infinity® Workstation** is the ultimate solution for today's classroom. Whether it is STEM, chemistry, biology, physics or a multitude of other disciplines, the Axis Infinity® is truly a "collaboratory" work space. The table is available in many different configurations and sizes to fit your needs.

sales@sheldonlabs.com



800.531.7604

AXIS INFINITY®

AXIS TABLE® MATRIX SHELDON LABORATORY SYSTEMS	MODEL #	REVIT FAMILY	TOP SIZE	HEIGHT CONTROL	SINK	RODS	SINK BASE	STUDENT CAPACITY
	26000	Axis Infinity	50" x 96" Epoxy	Motorized Adjustable 30-36" H	Infinity Sink®	Yes	Steel	4-6
	26200	Axis Infinity	44" x 84" Epoxy	Motorized Adjustable 30-36" H	Infinity Sink®	Yes	Steel	4
	26020	Axis Infinity No Sink	50" x 96" Epoxy	Motorized Adjustable 30-36" H	No	Yes	Steel	4-6
	26220	Axis Infinity No Sink	44" x 84" Epoxy	Motorized Adjustable 30-36" H	No	Yes	Steel	4
	26100	Axis Infinity	50" x 96" Epoxy	Motorized Adjustable 30-36" H	Infinity Sink®	Yes	Wood	4-6
	26300	Axis Infinity	44" x 84" Epoxy	Motorized Adjustable 30-36" H	Infinity Sink®	Yes	Wood	4
	26120	Axis Infinity No Sink	50" x 96" Epoxy	Motorized Adjustable 30-36" H, Wood Cover	No	Yes	No	4-6
	26320	Axis Infinity No Sink	44" x 84" Epoxy	Motorized Adjustable 30-36" H, Wood Cover	No	Yes	No	4
	26400	Axis Infinity Fixed Height	50" x 96" Epoxy	Fixed Ht. Adjust to 30, 34 or 36" H. Can retro. w/lift	#A25 Sink	No	Wood	4-6
	26500	Axis Infinity Fixed Height	44" x 84" Epoxy	Fixed Ht. Adjust to 30, 34 or 36" H. Can retro. w/lift	#A25 Sink	No	Wood	4
	26450	Axis-Synergy Sink-Fixed Height	50" x 96" Epoxy	Fixed Ht. Adjust to 30, 34 or 36" H. Can retro. w/lift	Infinity Sink®	Optional	Wood	4-6
	26550	Axis-Synergy Sink-Fixed Height	44" x 84" Epoxy	Fixed Ht. Adjust to 30, 34 or 36" H. Can retro. w/lift	Infinity Sink®	Optional	Wood	4
	26420	Axis-Synergy Fixed Height No Sink	50" x 96" Epoxy	Fixed Ht. Adjust to 30, 34 or 36" H. Can retro. w/lift	No	No	Wood	4-6
	26520	Axis-Synergy Fixed Height No Sink	44" x 84" Epoxy	Fixed Ht. Adjust to 30, 34 or 36" H. Can retro. w/lift	No	No	Wood	4
	26600	Axis Infinity- A25	50" x 96" Epoxy	Motorized Adjustable 30-36" H	#A25 Sink	No	Wood	4-6
	26700	Axis Infinity- A25	44" x 84" Epoxy	Motorized Adjustable 30-36" H	#A25 Sink	No	Wood	4
	26800 Physics	Axis Infinity- A25	50" x 96" Epoxy	Motorized Adj. 30-36" H Removable Sink Cover	#A25 Sink	No	Wood	4-6
	26900 Physics	Axis Infinity- A25	44" x 84" Epoxy	Motorized Adj. 30-36" H Removable Sink Cover	#A25 Sink	No	Wood	4
	Slim Line	Contact Sheldon Labs for more details on the SLIM LINE table options						
	26085 Mobile AXIS	Axis Infinity Mobile Sink	50" x 96" Epoxy	Motorized Adj. w/Casters 30-36" H. Includes: Sink w/Reservoir below	Removable 2-gal. Sink	Yes	Steel	4-6
	26275 Mobile AXIS	Axis Infinity Mobile Sink	44" x 84" Epoxy	Motorized Adj. w/Casters 30-36" H. Includes: Sink w/Reservoir below	Removable 2-gal. Sink	Yes	Steel	4
	26080 Mobile AXIS	Axis Infinity No Sink	50" x 96" Epoxy	Motorized Adjustable 30-36" H, w/casters	No	Yes	Steel	4-6
	26270 Mobile AXIS	Axis Infinity No Sink	44" x 84" Epoxy	Motorized Adjustable 30-36" H, w/casters	No	Yes	Steel	4
	21070	Harmony Docking Table	56"x28"x34" H Epoxy	Mobile Teaching Table with top detail to allow integral fit to docking at AXIS Infinity Table. Includes equipment shelf and leaf shelves on both sides (increases by table 24")				
	21075	Harmony Table	56"x28"x34" H Epoxy	Mobile Teaching Table with equipment shelf				
	Infinity Sink® Molded epoxy resin sink, 19 gals, 2 tiered with drain grooves, ant-splash feature, fits large dissection trays			#A25 Sink Molded epoxy resin sink- ADA 18" W x 15" L x 4.5" D			Privacy Panel #PPA0001 - 14" H x 19" D-Polymer Snaps to Axis Infinity Table to provide a stable vertical barrier that assures student privacy. #PPA0001K - 3 Panel Kit	
						Privacy Panel Storage Unit TM-426122 - Mobile unit that provides storage for 24 privacy panels and 14 Fiberglass Tote Trays #1700-0112. Outlet Options: A variety of electrical, ethernet & USB options are available		

sales@sheldonlabs.com



800.531.7604