



SYNERGY SINK

The **Synergy Sink** can be placed in the perimeter casework or in the Sheldon Axis Infinity® Table. With its curvilinear front, several students are able to collaborate simultaneously on their science project. The fixtures are mounted to have anti-splash characteristics directing the water flow to the sloped portion of the sink. This sink is versatile as it is useful.

Features

- 19 Gallons – fits large dissection trays
- Epoxy Resin
- Two-tiered sink with drain grooves
- Excellent for all disciplines of science
- ADA accessible is available

sales@sheldonlabs.com




800.531.7604

SYNERGY SINK



Synergy Sink (ADA) in perimeter casework



Synergy Sink in perimeter casework



Synergy Sink in Axis Infinity Table

Accessibility

Sink can be mounted at ADA height and incorporate compliant casework and metal knee panel. The front ADA section of the sink is 4.75" deep.

Fixtures

Unit comes standard with Unicast Cold Water/Gas Combination. Unimix Hot and Cold Water/Gas combination available at additional cost.

Sink Shape

Sink is two-tiered trapezoidal shape. Rear section of sink is 11.5" deep and 22.5" wide. Front section of sink is 4.75" deep and 28.375" wide.

Colors

Standard epoxy color is black. Limited variety of colors available at additional charge.

sales@sheldonlabs.com



800.531.7604