

A. **Sheldon 21075 Harmony Mobile Table:**
58" L. X 28" W. X 34" H.

1. The Harmony table allows mobile instruction within the classroom. The table is provided with a Shelresin top.
2. Four (4) large heavy-duty swivel casters (2-locking) support the table and large experiment storage shelf. The table has standard Sheldon thru-bolt corner design for complete stability.
3. The equipment storage shelf is constructed from 1" veneer core plywood mounted on leg stretchers for sound rigidity with 9" tall containment sides and back, and a front lip to keep experiment supplies stable. The table is available in Oak or Maple, as specified.
4. Width of the table allows transition from prep room storage to classroom easily through a 32" door. Rod sockets and rods are provided on the table to allow set-up prior to or during instruction.
5. Standard accessories to be included with each table:

One (1) 86380 Upright rod assembly.
Two (2) 86320 Rod bases.