

A. **Sheldon 31070 Potting Station:**
61" L. x 24" W. x 30" H. to work top.

TP 01.18.18

1. Work top is Phenolic Resin, 1" thick with 2" high curb at rear.
2. Hardwood upper shelf assembly to be 20-3/4" high x 8" deep and to have two (2) fixed shelves with rear curbs and vertical support ends.
3. One (1) B182922-100 base cabinet to have one adjustable shelf.
4. Knee space supported with leg assembly and fitted with two (2) 38802 soil carts with hardwood exterior and metal lined, water-tight interior. Cart to have hinged top with built-in stop.