

DP 08.12.20

A.

Sheldon P8201K Adjustable Height Table:

54"W x 24"D x 27-3/4"H to 39"H (42"H with included spacers).

1. The table provides a solution for ADA requirements that can adjust in increments from 27-3/4" to 39" high. It can also be adjusted for use as an incline surface for Physics. The solid hardwood apron frame is provided with a chemical resistant finish. The table has a weight capacity of 500 lbs.
2. The radius edged apron rails are 13/16" x 4-1/2" solid Oak and each securely attached with four (4) screws to curved, heavy-duty, steel corner braces. The interior face of aprons are provided with pocket holes for attaching frame to countertop with screws. Each adjustable height leg is attached to underside of the tabletop through a round, welded, heavy-duty steel plate with six (6) screws. Holes are pre-drilled for installation.
3. The unique design of the round silver powder coated steel Koyo legs allow for adjustment of the table height in less than 30 seconds without the use of tools.
4. The P8201K table is provided with following:
 - a. 1-1/4" Black Plastic Laminate Top.
 - b. Adjustable glides.
 - c. 3" spacers for additional height up to 42".
5. Available options as specified:
 - a. Casters (set of 4).
 - b. Maple apron frame.

Sheldon Labs

SheldonLabs.com
Sales@SheldonLabs.com
800.531.7604

Diversified Casework

DiversifiedCasework.com
Info@DiversifiedCasework.com
877.348.9663

LabsXpress

Info@LabsXpress.com