

DP 08.12.20

A. **Sheldon WB12-0V Four-Station Steel Workbench:**
64"W x 54"D x 33-1/4"H.

1. The 36"W x 42"D x 31"H locker base consists of a double-faced heavy gauge steel unit, welded and riveted throughout. The unit provides (12) vertical locker openings which measure 12"W x 21"D x 15"H. Vented doors have spring loaded hinges, padlock hasp, and knockout plug for cylinder locks. The base is provided with a steel angle bracket and two (2) gussets on each end to support top.
2. The 64" x 54" top is constructed of 2-1/4" thick maple.
3. The steel base is provided in gray color.
4. Available Options as specified:
 - a. WB12BL-0V with black steel base with no vises.
 - b. WB12-4V with four (4) #234708 woodworker's vises.
 - c. WB12BL-4V with black base and four (4) woodworker's vises.

Sheldon Labs

SheldonLabs.com
Sales@SheldonLabs.com
800.531.7604

Diversified Casework

DiversifiedCasework.com
Info@DiversifiedCasework.com
877.348.9663

LabsXpress

Info@LabsXpress.com