

DP 08.12.20

A. **Sheldon PT-60P Planning Table:**
60" W x 42" D x 30" H.

1. The PT Planning Table is designed to be used in most any environment from Art and General Planning to Sewing and Model construction. The base is constructed of two (2) tapered solid maple pedestal style panel legs with glides, and two (2) solid maple truss bars. Pedestal legs are supported with a horizontal solid maple stretcher rail. The top is uniquely attached to prevent the table from racking or rocking. The table is finished with a chemical resistant finish. Standard table height is 30".
2. The table is provided with:
A set of tote trays provided with tray runners and (4) totes.
Each tote measures 13-1/2"W x 19"D x 3-1/2"H.
1-1/4" almond plastic laminate top w/3mm PVC edge.
3. Table base is provided in Maple only.
4. Other Available Options as specified.

Sheldon Labs

SheldonLabs.com
Sales@SheldonLabs.com
800.531.7604

Diversified Casework

DiversifiedCasework.com
Info@DiversifiedCasework.com
877.348.9663

LabsXpress

Info@LabsXpress.com