

DP 08.12.20

A. **Sheldon P-7209K30N STEM Imprint Table:**  
54"W x 24"D x 30"H.

1. The Stem tables are built similar to plain apron tables, but with an easy-to-use imprinted top. The top is imprinted with graphs and measurement tools so students can measure and record shared data. The top has a compass, meter stick, yard stick, circular graph and data recording graph. The high-pressure laminate top works as a standard marker board, and can be marked up with crayons, grease pencil, or dry erase markers. Clean up is easy with standard markerboard cleaner or a damp cloth. The 54"W x 24"D top is perfect for up to four students to collaborate.
2. The table frames and legs are provided in solid hardwood. The radius edged apron rails are 13/16" x 4-1/2" and secured at each corner with 13 gauge, heavy-duty steel corner braces locked in grooves in rails, and attached with four (4) screws. The interior face of aprons are provided with pocket holes for attaching frame to countertop with screws. Legs are 2-1/4" square and radius edged, and secured to corner brackets with 5/16" hanger bolt, lock washer, and hex nut. Tables are provided with leg boots and levelers.

The table frame and legs are provided with a chemical resistant finish, and provided in Red Oak or Maple, as specified.

**Sheldon Labs**

SheldonLabs.com  
Sales@SheldonLabs.com  
800.531.7604

**Diversified Casework**

DiversifiedCasework.com  
Info@DiversifiedCasework.com  
877.348.9663

**LabsXpress**

Info@LabsXpress.com