

DP 04.02.20

A. **Sheldon AFT-483630-B-LO9-09-C A-Frame Table:**  
48" W. X 36" D. X 30" H. (Nominal height depending on top thickness).

1. The A-Frame Table is a versatile mobile table provided with a stylish heavy gauge steel frame with full steel tube legs, additional welded cross braces, formed top leg brackets, and formed C-channel stringer brackets. The steel frame and legs have a baked powder coat finish.
2. Each leg frame assembly has two (2) 14-ga., 2-1/2" x 2-1/2" steel tube legs with welded bottom plate and threaded nut, one (1) 14-ga., 2" x 2" steel tube leg cross brace welded to legs, and one (1) 11-ga., 2-3/4" x 2" L-shaped top leg bracket welded to legs.
3. Leg frame assemblies are secured to two (2) 14-ga., 2-1/2" x 2" formed steel C-channel stringers attached each end with two (2) bolts thru the stringer end plate, leg, and top leg bracket.
4. This table frame is provided with:  
Black (B) finish.  
1-1/4" Black Laminate Top (L09) with Black Edge Band (09).  
4" Industrial Casters (C).
5. Available Options as specified:
  - a. Width: 48", 60", 72".
  - b. Depth: 36", 42", 48". (Frame depth is 34").
  - c. Height: 30", 36".
  - d. Caster and Glide Options:  
4"Industrial Caster (C), 4"Designer Caster (D),  
Leveling Glide (G).
  - e. Frame colors: Black (B), Silver (S), White (W),  
Orange (O).
  - f. Top options: 1" Epoxy (E), 3/4" Phenolic (P), 1-  
1/2" Shop Top (S), 1-3/4" Butcher Block -  
Walnut (BBW) and White (BBN).
  - g. 1-1/4" Laminate Tops - Black (L09), Wrought  
Iron (LL07),

#### Sheldon Labs

SheldonLabs.com  
Sales@SheldonLabs.com  
800.531.7604

#### Diversified Casework

DiversifiedCasework.com  
Info@DiversifiedCasework.com  
877.348.9663

#### LabsXpress

Info@LabsXpress.com



- h. Almond (L04), and White (L03).  
Edge Band Options for Laminate Tops - Black (09), Almond (04), Wrought Iron (L07), White (03), Papaya (C98), Island (A40), Ocean (A41), Lapis Blue (14), Hollyberry (60), and Summer Vacation (D13).

**Sheldon Labs**

SheldonLabs.com  
Sales@SheldonLabs.com  
800.531.7604

**Diversified Casework**

DiversifiedCasework.com  
Info@DiversifiedCasework.com  
877.348.9663

**LabsXpress**

Info@LabsXpress.com