



HOT WATER MOBILE STATION
HSW-3624K

sheldonlabs™

PRODUCT

KEY FEATURES

Hot Water Mobile Station

- Oak veneer, locking double door cabinet
- Black plastic laminate top
- Water heater and pump provides hot or cold water
- 5 gallon fresh water and 6 gallon waste water tanks are NSF, UL and FDA approved
- Stainless steel sink
- Hot and cold water mixing faucet
- Duplex outlet with included extension cord
- For optimal performance start with room temperature water
- Water heats up to 10 to 15 degrees warmer than the fresh tank temperature
- Rolls easily on 4" locking, swivel casters

ITEM#	PRODUCT	SIZE
HSW-3624	Hot Water Mobile Station	36"W x 24"D x 36"H