

TP 01.29.20

A. **Sheldon IMW-92K-36H-L04-04 Intermix Mobile Workstation:**

30" W. X 30" D. X 36" H. Lazy Susan Cabinet with 60" diagonal overall top.

1. The IMW-92 Intermix Workstation is a mobile storage unit consisting of a dowel constructed Oak plywood open cabinet provided with a chemical resistant UV finish, and mounted on four (4) locking swivel casters. The four-sided unit contains a compartment with a Lazy Susan unit having two (2) independent turntables accessible from all four (4) sides. The cabinet has a 3/4" top panel, bottom panel, and four (4) side panels, along with an additional 1" applied bottom panel to provide support for attachment of the four (4) casters.
2. This mobile unit has options of Lazy Susan - (92), Oak Cabinet - (K), 36" High - (36H), 1-1/4" Almond Laminate Top - (L04), and Almond edge band - (04).
3. Available Options as specified:
 - a. Height: 36" High - (36H) or 30" High - (30H).
 - b. Cabinet: Oak - (K) or Maple - (M).
 - c. Countertops: 3/4" Phenolic Resin - (P), 1-1/2" Shop Top - (S), 1-1/4" Laminates: Black - (L09), Wrought Iron - (LL07), Almond - (L04), and White - (L03).
 - d. Edge band choices: Black - (09), Wrought Iron - (L07), Almond - (04), White - (03), Papaya - (C98), Island - (A40), Ocean - (A41), Lapis Blue - (14), Hollyberry - (60), Summer Vacation - (D13).

Sheldon Labs

SheldonLabs.com
Sales@SheldonLabs.com
800.531.7604

Diversified Casework

DiversifiedCasework.com
Info@DiversifiedCasework.com
877.348.9663

LabsXpress

Info@LabsXpress.com