

TP 01.29.20

A. **Sheldon IMW-03K-36H-L04-04 Intermix Mobile Workstation:**

30" W. X 30" D. X 36" H. Storage Cabinet with 60" diagonal overall top.

1. The IMW-03 Intermix Workstation is a mobile storage unit consisting of a dowel constructed Oak plywood cabinet provided with a chemical resistant UV finish, and mounted on four (4) locking swivel casters. The double-sided unit has double doors with a lock. It has a storage compartment on front and rear sides with a 3/4" center divider panel separating the two (2) storage areas. The cabinet has a 3/4" top panel, bottom panel, and end panels, along with an additional 1" applied bottom panel to provide support for attachment of the four (4) casters. Each storage compartment has a 3/4" x 12-1/2" deep adjustable shelf.
2. This mobile unit has options of Doors - (03), Oak Cabinet - (K), 36" High - (36H), 1-1/4" Almond Laminate Top - (L04), and Almond edge band - (04).
3. Available Options as specified:
  - a. Height: 36" High - (36H) or 30" High - (30H).
  - b. Cabinet: Oak - (K) or Maple - (M).
  - c. Countertops: 3/4" Phenolic Resin - (P), 1-1/2" Shop Top - (S), 1-1/4" Laminates: Black - (L09), Wrought Iron - (LL07), Almond - (L04), and White - (L03).
  - d. Edge band choices: Black - (09), Wrought Iron - (L07), Almond - (04), White - (03), Papaya - (C98), Island - (A40), Ocean - (A41), Lapis Blue - (14), Hollyberry - (60), Summer Vacation - (D13).

**Sheldon Labs**

SheldonLabs.com  
Sales@SheldonLabs.com  
800.531.7604

**Diversified Casework**

DiversifiedCasework.com  
Info@DiversifiedCasework.com  
877.348.9663

**LabsXpress**

Info@LabsXpress.com