

A. **Sheldon IMW-99K-36H-L04-04 Intermix Mobile Workstation:**

30" W. X 36" D. X 36" H. Tote Tray Cabinet with 60" diagonal overall top.

1. The IMW-99 Intermix Workstation is a mobile storage unit consisting of a dowel constructed Oak plywood open cabinet provided with a chemical resistant UV finish, and mounted on four (4) locking swivel casters. The double-sided unit contains two (2) compartments on both front and rear sides, and a 3/4" center divider panel separates the front and rear compartments. The cabinet has a 3/4" top panel, bottom panel, and end panels, along with an additional 1" applied bottom panel to provide support for attachment of the four (4) casters.
2. The front and rear compartments each have a 3/4" vertical divider panel that separates the compartment into two (2) sections. One section has a fixed shelf for storage. The other section has a book compartment above and a space below that holds five (5) Certwood plastic tote trays with pull-out glide and tilt runners. This space includes four (4) small tote trays 12-1/4" W. x 16-3/4" D. x 3" H., and one (1) larger tote tray that is 12-1/4" W. x 16-3/4" D. x 6" high.
 - a. The Certwood Tote trays are shatterproof, clear ABS plastic designed with a dual labeling system. The front side has a removable translucent ticket window that provides the option to utilize an internal label or external writable surface. The rear side has a 2" high x 8-1/4" wide "E-Z Peel" area with stipples no less than .004 above the face of the tray, providing a surface to easily apply and remove larger sticker labels. Trays have a 3/16" internal radius for ease of cleaning.
 - b. The Glide & Tilt Tray Runner System is also constructed of clear ABS. The runners hold trays in place and allow them to be opened in an angled position for easy access to contents.

Sheldon Labs

SheldonLabs.com
Sales@SheldonLabs.com
800.531.7604

Diversified Casework

DiversifiedCasework.com
Info@DiversifiedCasework.com
877.348.9663

LabsXpress

Info@LabsXpress.com

Trays are also reversible for use of the front or back dual labeling system.

- c. The triple module tote tray runners are designed with integral tab inserts that inset into the pre-drilled vertical cabinet end panel and divider panel, and are attached front and back with screws.
3. This unit has options of Totes - (99), Oak Cabinet - (K), 36" high - (36H), 1-1/4" Almond Laminate Top - (L04), and Almond edge band - (04).
 4. Available Options as specified:
 - a. Height: 36" High - (36H) or 30" High - (30H). (30" high unit has (2) small totes & (1) larger tote).
 - b. Cabinet: Oak - (K) or Maple - (M).
 - c. Tops: 3/4" Phenolic - (P), 1-1/2" Shop Top - (S), 1-1/4" Laminates: Black - (L09), Wrought Iron - (LL07), Almond - (L04), White -(L03).
 - d. Edge band: (09), (L07), (04), (03), (C98), (A40), (A41), (14), (60), and (D13).

Sheldon Labs

SheldonLabs.com
Sales@SheldonLabs.com
800.531.7604

Diversified Casework

DiversifiedCasework.com
Info@DiversifiedCasework.com
877.348.9663

LabsXpress

Info@LabsXpress.com