

4-Door,  
Lazy Susan,  
and Tote  
Options



**INTERMIX WORKSTATION**  
IMW-99K-36H-LW-R



**INTERMIX WORKSTATION**  
IMW-99M-30H-LC-A

sheldonlabs™

**PRODUCT**

**KEY FEATURES**

**Intermix Workstation**

- 2 table heights, 30" & 36"
- Unique top shape fosters collaboration and encourages productivity
- Available in 3 base configurations: Lazy-Susan, 4-Door, and Open Tote Cabinets
- Lazy-Susan Cabinet with two independent turntables
- 4-Door Cabinet with an adjustable shelf and door locks
- Open Tote Cabinet (36"H): 1 double and 4 single totes per side
- Open Tote Cabinet (30"H): 1 double and 2 single totes per side
- Single totes (12-1/4"W x 16-3/4"D x 3"H)
- Double totes (12-1/4"W x 16-3/4"D x 6"H)
- Solid oak or maple veneers finished with a chemical resistant, earth-friendly UV finish
- 1-1/4" high pressure laminate tops with durable 3mm PVC edge band
- 4 locking swivel casters for ease of mobility
- Durable, plastic laminate top with 4 color options
- 10 edge-banding options for laminate tops
- Overall top size: 60"
- Assembly required

**SIZE OPTIONS**

- Width: 60"
- Depth: 60"
- Height: 30", 36"

**WOOD BASE OPTIONS**

- K: Oak
- M: Maple

**CABINET OPTIONS**

- 03: 4-Door Cabinet
- 92: Lazy Susan
- 99: Open Tote

**TOP OPTIONS**

- LA: 1-1/4" Almond Laminate\*
- LB: 1-1/4" Black Laminate\*
- LC: 1-1/4" Charcoal Laminate\*
- LW: 1-1/4" White Laminate\*

**\* EDGE BANDING OPTIONS FOR LAMINATE TOPS**

(Must specify edge banding if laminate top is a different color)

- A: Almond
- B: Black
- C: Charcoal
- L: Lime
- P: Peacock
- Q: Aqua
- R: Red
- T: Carrot
- W: White
- Y: Royal

**BASE CABINET  
OPTIONS**



**Lazy-Susan**  
IMW-92M-30H-LA




**4-Door Cabinet**  
IMW-03M-30H-LA




**Open Tote Cabinet**  
IMW-99M-30H-LC-A

**IMW-99K-30H-LA-Q**




**IMW-99K-30H-LA-Q**  
30"H, Open Tote Cabinet Base, Oak,  
Almond Top, Aqua Edge Banding

**PLASTIC LAMINATE  
TOP OPTIONS**



**PLASTIC LAMINATE  
EDGE BANDING  
OPTIONS**



MODEL	WOOD	HEIGHT	BLACK LAMINATE	ALMOND LAMINATE	WHITE LAMINATE	CHARCOAL LAMINATE
4 Door	Oak	30"	IMW-03K-30H-LB	IMW-03K-30H-LA	IMW-03K-30H-LW	IMW-03K-30H-LC
4 Door	Oak	36"	IMW-03K-36H-LB	IMW-03K-36H-LA	IMW-03K-36H-LW	IMW-03K-36H-LC
4 Door	Maple	30"	IMW-03M-30H-LB	IMW-03M-30H-LA	IMW-03M-30H-LW	IMW-03M-30H-LC
4 Door	Maple	36"	IMW-03M-36H-LB	IMW-03M-36H-LA	IMW-03M-36H-LW	IMW-03M-36H-LC
Lazy Susan	Oak	30"	IMW-92K-30H-LB	IMW-92K-30H-LA	IMW-92K-30H-LW	IMW-92K-30H-LC
Lazy Susan	Oak	36"	IMW-92K-36H-LB	IMW-92K-36H-LA	IMW-92K-36H-LW	IMW-92K-36H-LC
Lazy Susan	Maple	30"	IMW-92M-30H-LB	IMW-92M-30H-LA	IMW-92M-30H-LW	IMW-92M-30H-LC
Lazy Susan	Maple	36"	IMW-92M-36H-LB	IMW-92M-36H-LA	IMW-92M-36H-LW	IMW-92M-36H-LC
Open Tote	Oak	30"	IMW-99K-30H-LB	IMW-99K-30H-LA	IMW-99K-30H-LW	IMW-99K-30H-LC
Open Tote	Oak	36"	IMW-99K-36H-LB	IMW-99K-36H-LA	IMW-99K-36H-LW	IMW-99K-36H-LC
Open Tote	Maple	30"	IMW-99M-30H-LB	IMW-99M-30H-LA	IMW-99M-30H-LW	IMW-99M-30H-LC
Open Tote	Maple	36"	IMW-99M-36H-LB	IMW-99M-36H-LA	IMW-99M-36H-LW	IMW-99M-36H-LC