

## DP 10.29.20

A. **Sheldon AFT-603630-S-LB-A-C A-Frame Table:**  
60" W. X 36" D. X 30" H. (Nominal height depending on top thickness).

1. The A-Frame Table is a versatile mobile table provided with a stylish heavy gauge steel frame with full steel tube legs, additional welded cross braces, formed top leg brackets, and formed C-channel stringer brackets. The steel frame and legs have a baked powder coat finish.
2. Each leg frame assembly has two (2) 14-ga. 2-1/2" x 2-1/2" steel tube legs with welded bottom plate and threaded nut, one (1) 14-ga. 2" x 2" steel tube leg cross brace welded to legs, and one (1) 11-ga. 2-3/4" x 2" L-shaped top leg bracket welded to legs.
3. Leg frame assemblies are secured to two (2) 14-ga. 2-1/2" x 2" formed steel C-channel stringers attached each end with two (2) bolts through the stringer end plate, leg, and top leg bracket.
4. This AFT A-Frame Table has options of:
  - a. Width: 60", Depth: 36", Height: 30".
  - b. Frame color - S: Silver.
  - c. Top - LB: 1-1/4" Black Laminate.
  - d. 3mm PVC Edge Band - A: Almond.
  - e. C: 4" Industrial Casters.
5. See **sheldonlabs** Sales Sheet for available options that may be Specified.

Frame Sizes - width, depth, and height.  
Frame colors.  
Casters & Leveling Glide.  
Tops materials.  
Edge Banding colors for Laminate tops.

### Sheldon Labs

SheldonLabs.com  
Sales@SheldonLabs.com  
800.531.7604

### Diversified Casework

DiversifiedCasework.com  
Info@DiversifiedCasework.com  
877.348.9663

### LabsXpress

Info@LabsXpress.com