



Height Adjustable
28" to 45"H



Fold-Away Hand Crank

sheldonlabs™

PRODUCT

KEY FEATURES

Hi-Lo Bench

- Designed for ergonomic and ADA needs
- Work surface adjusts for a comfortable and efficient use
- Height adjusts from 28" to 45"H
- Turning the hand crank 2 times moves the work surface 1"
- The fold-away hand crank stores conveniently under the table
- Add optional casters for mobility
- Top options: 1-3/4" solid maple; 1" epoxy resin; 1-1/4" plastic laminate top in either black or almond
- Frame color: Black

ITEM#	TOP	SIZE
HLB-6030M	Maple Butcher Block	60"W x 30"D x 28"-45"H
HLB-6030PA	Plastic Laminate - Almond	60"W x 30"D x 28"-45"H
HLB-6030PB	Plastic Laminate - Black	60"W x 30"D x 28"-45"H
HLB-7224ER	Epoxy Resin	72"W x 24"D x 28"-45"H
HLB-7224M	Maple Butcher Block	72"W x 24"D x 28"-45"H
HLB-7224PA	Plastic Laminate - Almond	72"W x 24"D x 28"-45"H
HLB-7224PB	Plastic Laminate - Black	72"W x 24"D x 28"-45"H
HLB-7230ER	Epoxy Resin	72"W x 30"D x 28"-45"H
HLB-7230M	Maple Butcher Block	72"W x 30"D x 28"-45"H
HLB-7230PA	Plastic Laminate - Almond	72"W x 30"D x 28"-45"H
HLB-7230PB	Plastic Laminate - Black	72"W x 30"D x 28"-45"H
253941	Casters (Set of 4)	---