

A. **Sheldon HWS-3624K Hot Water Mobile Station:**
36"W x 24"D x 36"H

1. The Hot Water Mobile Station includes multiple features.
 - a. Oak veneer, locking double door cabinet.
 - b. 1-1/4" Black plastic laminate top.
 - c. Water heater and pump provides hot or cold water.
 - d. Outside pump switch and inside water heater switch.
 - e. 5 gallon freshwater and 6 gallon wastewater tanks (NSF, UL, and FDA approved).
 - f. Stainless steel sink.
 - g. Hot and cold water mixing faucet.
 - h. Duplex outlet with included extension cord.
 - i. For optimal performance start with room temperature water.
 - j. Water heats up to 10 to 15 degrees warmer than the fresh tank temperature.
 - k. Rolls easily on 4" locking, swivel casters.
 - l. Lockable doors.

2. Oak Veneer Cabinet.
 - a. 34"W x 29.50" high x 22"D. Without laminate top and casters.
 - b. 3/4" doweled constructed, oak veneer end, back, bottom and door panels, banded on exposed edges.

3. See **sheldonlabs** Sales Sheet for description and features.