

Contour  
Seat  
New!



BLUE - STANDING HEIGHT  
WITH OPTIONAL CASTERS  
SE-TR1M



BLACK - DESK HEIGHT  
SE-TR2D

sheldonlabs™

PRODUCT

KEY FEATURES

## Tractor Stools

- Lab compliant clean room seating
- Self-skinning polyurethane that is both chemical and puncture resistant for exceptional durability
- Contoured tractor style for supreme comfort
- Fingertip pneumatic height adjustment
- Height adjustable foot ring on standing models
- 25" glass fiber reinforced nylon base
- Non-marring glides standard, self-braking casters optional
- Rated for 300# weight capacity
- Seat size 17"W x 14"D
- Glass reinforced nylon base
- 15 year warranty
- Minimal assembly

COLOR OPTIONS

- Black
- Blue

HEIGHT OPTIONS

- Desk height range 16" – 21"
- Standing height range 23" – 33"



SELF-BRAKING CASTERS (SE-C)

ITEM #	PRODUCT	COLOR	ADJUSTABLE HEIGHT	FOOT RING	SEAT SIZE
SE-TR2D	Tractor Stool - Desk	Black	16" - 21"H	No	17"W x 14"D
SE-TR2M	Tractor Stool - Standing	Black	23" - 33"H	Yes	17"W x 14"D
SE-TR1D	Tractor Stool - Desk	Blue	16" - 21"H	No	17"W x 14"D
SE-TR1M	Tractor Stool - Standing	Blue	23" - 33"H	Yes	17"W x 14"D
SE-C	Castors (Qty: 5)	Black	--	--	--