




Hydroponics Growing Center
HGC-52



sheldonlabs™

PRODUCT

KEY FEATURES

Hydroponics Growing Center

Self-contained ebb and flow hydroponics growing system. The perfect way to learn how plants can grow without soil and helps teach the basic concepts for the plants needs, nutrition, food production, recycling, agricultural technology and soil-less growing environments. This mobile hydroponics center includes:

- Grow lights and water pump
- Adjustable timers enable you to grow your plants with the perfect amount of water and sunlight.
- Constructed of solid maple and maple veneers
- Includes containment plant tray (8"H x 27-1/2"D x 42"W)
- 25 gallon nutrient reservoir with lid
- Locked cabinet for timers and air pump
- 150 watt grow light.
- Submersible pump, air pump, air stone, (2) timers, hoses and ten 6" net pots.
- Moves easily on 4" locking, leveling casters.
- Overall dimensions: 52"W x 28"D x 78"H

ITEM #	PRODUCT	SIZE
HGC-52	Hydroponics Growing Center	52"W x 28"D x 78"H