

- A. **Sheldon 260004 Cold Water Sink & Metered Faucet Mobile Unit:**
24"W x 24"D x 38.81"H
1. The Mobile Handwash Sink is made of heavy duty stainless steel construction and includes multiple features.
 - a. 4-gauge stainless steel top.
 - b. 18-gauge stainless steel body.
 - c. Four locking 5" casters that allow for easy mobility (2-swivel and 2-rigid).
 - d. Lockable access door.
 - e. Bowl size: 10" x 14" x 5".
 - f. Water-saving metering faucet.
 - g. 5-gallon removable freshwater tank.
 - h. No disassembly required to refill.
 - i. 7-gallon removable wastewater tank.
 - j. Drain hose attached for easy drain.
 - k. Soap and towel dispenser.
 - l. Top fill soap dispenser for easy refill.
 - m. 115 volt 3.3 GPM, 45 PSI water diaphragm pressure pump (CSA listed).
 - n. 115 volt water heater – 2.5 gallons (ETL listed), in hot water units only.
 - o. 15 ft. GFCI power cord with standard wall plug (UL listed).
 - p. Hose connections provided for water supply and drain connections.
 2. See **sheldonlabs** Sales Sheet for available options that may be specified for Hands free electronic sensor faucet, or for hot water only units.