

A. **Sheldon HLB-6030M Hi-Lo Bench:**
60"W x 30"D with adjustable height.

1. The Hi-Lo Bench is designed for ergonomic and ADA needs, and allows the work surface to be adjusted for a comfortable and efficient positioning of the arms and hands.
2. The black steel table frame is provided with leveling glides, and can be adjusted to a height of 26-1/2" to 39-3/4" from floor to underside of optional work surfaces. Add work surface thickness for total height adjustment including glides. Add additional 1-3/4" to table height if optional casters for mobility are included.
3. The fold-away hand crank stores conveniently under the table. Turning the hand crank (2) times moves the work surface 1 inch.
4. The HLB-6030M table is provided with the following:
 - a. 1-3/4" Maple Butcher Block top. 60"W x 30"D.
 - b. Adjustable glides.
 - c. Height adjustment with glides and Maple top is 41-1/2".
5. See **sheldonlabs** Sales Sheet for available sizes, work surface options, and casters (set of 4) that may be specified.